

Product Datasheet

Mouse PRSS27 ELISA Kit (orb385396)

Description Sandwich ELISA for quantitative detection of mouse Pancreasin concentrations.

Reactivity Mouse

Range 125 pg/ml - 8000 pg/ml

Tested Applications ELISA

Target PRSS27

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum

Uniprot ID [Q8BJR6](#)

Sensitivity 10 pg/ml

Entrez [213171](#)

Expiration Date Please enquire.



Representative
Standard
Curve