

## Product Datasheet

### Mouse TGFBI ELISA Kit (orb385403)

**Description** Sandwich ELISA for quantitative detection of mouse TGFBI concentrations.

**Reactivity** Mouse

**Range** 156 pg/ml - 10000 pg/ml

**Tested Applications** ELISA

**Target** TGFBI

**Storage** Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

**Note** For research use only

**Application notes** Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


**Sample Types** Cell Culture Supernates, Serum, Plasma

**Uniprot ID** **P82198**

**Sensitivity** 5 pg/ml

**Entrez** **21810**

**Expiration Date** Please enquire.



Representative  
Standard  
Curve