

Product Datasheet

Mouse TNFRSF25 ELISA Kit (orb385405)

Description Sandwich ELISA for quantitative detection of mouse DR3 concentrations.

Reactivity Mouse

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target TNFRSF25

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **B1AWN9**

Sensitivity 15 pg/ml

Entrez **85030**

Expiration Date Please enquire.



Representative
Standard
Curve