

Product Datasheet

Mouse TREML1 ELISA Kit (orb385406)

Description Sandwich ELISA for quantitative detection of mouse TLT1 concentrations.

Reactivity Mouse

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target TREML1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q8K558**

Sensitivity 10 pg/ml

Entrez **71326**

Expiration Date Please enquire.



Representative
Standard
Curve