


Product Datasheet

Mouse WIF1 ELISA Kit (orb385408)

Description	Sandwich ELISA for quantitative detection of mouse WIF-1 concentrations.
Reactivity	Mouse
Range	23.4 pg/ml - 1500 pg/ml
Tested Applications	ELISA
Target	WIF1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9WUA1
Sensitivity	10 pg/ml
Entrez	24117
Expiration Date	Please enquire.



Representative
Standard
Curve