



Product Datasheet

MOGAT1 antibody (orb385431)


Description	MOGAT1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MOGAT1. Please contact us for the exact immunogen sequence. The peptide is available as orb13062.
Target	MOGAT1. This antibody has specificity with MOGAT1.
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	39 kDa
Uniprot ID	D4A6J3 , Q96PD6 , Q91ZV4
NCBI	058165 , 477513
Entrez	116255
Dilution Range	IHC-P: 1:100-800, IF/ICC: 1:100-800, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



IHC-P image of mouse liver tissue using ...



IHC-P staining of human lung cancer tiss...



IHC-P staining of human lung cancer tiss...