

Product Datasheet

ASL1A antibody (orb385432)

Description

ASL1A antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human ASL1A. Please contact us for the exact immunogen sequence. The peptide is available as orb13063.

Target

ASL1A. This antibody has specificity with ASL1A.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

25 kDa

Uniprot ID

[Q9Y294](#), [Q9CQE6](#)

NCBI

[014034](#), [054753](#)

Entrez


[25842](#)

Dilution Range


IHC-P: 1:50-300, WB: 1:200-2000

Expiration Date

12 months from date of receipt.



Western blot analysis of mouse brain (la...



Immunohistochemical staining of paraffin...