

Product Datasheet

Nesprin 1 antibody (orb385439)

Description

Nesprin 1 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Nesprin 1. Please contact us for the exact immunogen sequence. The peptide is available as orb13070.

Target

Nesprin 1. This antibody has specificity with Nesprin 1.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

112 kDa

Uniprot ID

[Q8VHJ9](#), [Q8NF91](#), [Q6ZWR6](#)

NCBI

[149062](#), [033071](#)

Entrez

[23345](#)


Dilution Range

IHC-P: 1:100-600, IF/ICC: 1:100-600, WB: 1:200-2000


Expiration Date

12 months from date of receipt.


IF IF IF



IF image of rat pancreas tissue using Ne...



IHC-P image of rat pancreas tissue using...



Immunofluorescence analysis of guinea pi...