

Product Datasheet

CD11A antibody (orb385445)

Description

CD11A antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD11A. Please contact us for the exact immunogen sequence. The peptide is available as orb13080.

Target

CD11A. This antibody has specificity with CD11A.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

130kDa

Uniprot ID

P24063, **F1LP44**, **P20701**

NCBI

001107852, **001114380**

Entrez


3683

Dilution Range


IHC-P: 1:100-1:600, WB: 1:100-1:800

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat lymph...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat lymph...