

Product Datasheet

CD169 Sialoadhesin antibody (orb385451)

Description

CD169 Sialoadhesin antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human CD169 Sialoadhesin. Please contact us for the exact immunogen sequence. The peptide is available as orb13102.

Target

CD169 Sialoadhesin. This antibody has specificity with CD169 Sialoadhesin.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

IHC-P staining of guinea pig ovary tissue...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P image of guinea pig bladder tissue...

MW

183 kDa

Uniprot ID
[Q62230](#), [Q9BZZ2](#)
NCBI
[075556](#), [023068](#)
Entrez
[6614](#)
Dilution Range

IHC-P: 1:100-500, IF/ICC: 1:100-500, WB: 1:200-1000

Expiration Date

12 months from date of receipt.