



Product Datasheet

ATP1A2 antibody (orb385453)


Description	ATP1A2 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ATP1A2. Please contact us for the exact immunogen sequence. The peptide is available as orb13104.
Target	ATP1A2. This antibody has specificity with ATP1A2.
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	112 kDa
Uniprot ID	Q6PIE5, P06686, P50993
NCBI	000702, 000693
Entrez	477
Dilution Range	IHC-P: 1:100-500, IF/ICC: 1:100-500, WB: 1:200-1500
Expiration Date	12 months from date of receipt.



IHC-P staining of guinea pig brain tissue...



IHC-P staining of mouse heart tissue using...



IHC-P image of rat prostate tissue using...