


Product Datasheet

NR1I2 antibody (orb385455)

Description

NR1I2 antibody



Species/Host


Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated



Tested

ICC, IF, IHC-P, WB

Applications

Immunohistochemical staining of paraffin...

Immunogen


KLH conjugated synthetic peptide derived from human NR1I2. Please contact us for the exact immunogen sequence. The peptide is available as orb13108.

Target

NR1I2. This antibody has specificity with NR1I2.

Form/Appearance

10 mM PBS, 0.02% sodium azide



Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Immunohistochemical staining of paraffin...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note


For research use only

Isotype

IgG

Clonality

Polyclonal



Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P image of human breast cancer tissue...

MW

50 kDa

Uniprot ID

[O75469](#), [O54915](#), [Q9R1A7](#)

NCBI

[003880](#), [003889](#)

Entrez

[8856](#)

Dilution Range

IHC-P: 1:100-600, IF/ICC: 1:100-600, WB: 1:200-1500

Expiration Date

12 months from date of receipt.