

Product Datasheet

Tyrosine Hydroxylase antibody (orb385456)

Description

Tyrosine Hydroxylase antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Tyrosine Hydroxylase. Please contact us for the exact immunogen sequence. The peptide is available as orb13111.

Target

Tyrosine Hydroxylase. This antibody has specificity with Tyrosine Hydroxylase.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

59 kDa

Uniprot ID

[P24529](#), [P07101](#), [P04177](#)

NCBI

[000351](#), [000360](#)

Entrez


[7054](#)

Dilution Range


IHC-P: 1:100-800, IF/ICC: 1:100-800, WB: 1:200-1500

Expiration Date


12 months from date of receipt.



IHC-P staining of guinea pig brain tissue...



Immunohistochemical staining of guinea p...



Immunohistochemical staining of mouse ad...