

Product Datasheet

Autotaxin antibody (orb385459)

Description

Autotaxin antibody

Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Autotaxin. Please contact us for the exact immunogen sequence. The peptide is available as orb394010.
Target	Autotaxin. This antibody has specificity with Autotaxin.
Form/Appearance	10 mM PBS, 0.02% sodium azide
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	25 kDa
Uniprot ID	Q9R1E6 , Q13822 , Q64610
NCBI	001040092 , 001035181
Entrez	5168
Dilution Range	IHC-P: 1:100-500, IF/ICC: 1:100-500, WB: 1:200-1500
Expiration Date	12 months from date of receipt.