


Product Datasheet

ACO2 antibody (orb385466)

Description	Goat polyclonal antibody to ACO2.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Aconitase 2 (aa541-555), Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in Human Heart lysates, in Human, Mouse and Rat Adipose lysates and in Mouse, Rat and Pig Skeletal Muscle lysates (calculated MW of 85.4kDa according to NP_001089.1)
Protein Sequence	CQDTYQHPPKDSSGQH.
Clonality	Polyclonal
Entrez	50
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western Blot of Human Heart lysate using...