


Product Datasheet

AP2A1 antibody (orb385471)

Description	Goat polyclonal antibody to AP2A1.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AP2A1 (aa706-727), Biotinlyated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 110kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 108kDa according to NP_055018.2)
Protein Sequence	CELEPPAPESPMALLADPAPAAD.
Clonality	Polyclonal
Entrez	160
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western Blot
of Human
Frontal Cortex
lys...