


Product Datasheet

BMP7 antibody (orb385477)

Description	Goat polyclonal antibody to BMP7.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BMP7 (aa120-131)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50-55kDa band observed in Pig Kidney lysates (calculated MW of 49.2kDa according to XP_005673101.1)
Protein Sequence	SLQDSHFLTDAD
Clonality	Polyclonal
MW	49.3
Entrez	655
Dilution Range	ELISA: 1:2000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western Blot
of Pig Kidney
lysate using

...