

Product Datasheet

KMO antibody (orb385480)

Description	Goat polyclonal antibody to KMO.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	kynurenine 3-monooxygenase (aa29-42)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human Placenta lysates (calculated MW of 55.8kDa according to NP_003670.2)
Protein Sequence	KRNFQIDVYEARED
Clonality	Polyclonal
MW	55.8
Entrez	8564
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western Blot
of Human
Placenta
lysate us...