


Product Datasheet

LCAT antibody (orb385481)

Description	Goat polyclonal antibody to LCAT.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	LCAT (aa366-378)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human Testis lysates (calculated MW of 49.6kDa according to NP_000220.1)
Protein Sequence	EDGDDTVATRSTE
Clonality	Polyclonal
MW	49.6
Entrez	3931
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western Blot
of Human
Testes lysate
usin...