

Product Datasheet

Interleukin-1 receptor-like 1 antibody (orb385505)

Description	Rabbit polyclonal antibody to Interleukin-1 receptor-like 1
Species/Host	Rabbit
Conjugation	Unconjugated
Concentration	See vial for concentration
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Isotype	Ig
Uniprot ID	Q01638
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.