

Product Datasheet

S-2238 (Thrombin Substrate) (orb385603)

Description	S-2238 is specific to thrombin, is a short peptide covalently bound to pNA (4-nitroaniline).
Note	For research use only
Formula	H-D-Phe-Pip-Arg-pNA·2HCl
MW	625.6
CAS Number	[115388-96-0]
SMILES	<chem>C1CCN(C(C1)C(=O)N[C@@H](CCCN=C(N)N)C(=O)NC2=CC=C(C=C2)[N+](=O)[O-])C(=O)[C@@H](CC3=CC=CC=C3)N</chem>
Expiration Date	12 months from date of receipt.