

Product Datasheet

HCG-beta Antibody (orb385662)

Description

This antibody reacts with a protein of 22kDa,

cross r...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-Fr

Applications

Immunogen

Purified HCG-beta protein was used as the immunogen for this antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the HCG-beta antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the HCG-beta antibody to be titrated up or down for optimal performance.1. Clone HCGb/211 (recommended capture) will pair with clones **HCGb/459 and HCGb/54 (recommended detect) by ELISA.**

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This HCG-beta antibody is available for research use only.

Entrez

1082

Dilution Range

Immunohistochemistry (Frozen): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.