

Product Datasheet

CD19 Antibody / B-lymphocyte marker (orb385675)

Description

CD19 is a transmembrane glycoprotein

immunoglobulin-like domains. CD...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

Human CD19 recombinant protein was used as immunogen for this CD19 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

CD19 antibody with azide can be stored at 2-8°C. The azide-free format should be aliquoted and stored at -20°C or colder.

Note

For research use only

Application notes

Variations in protocols, secondaries and substrates may require the CD19 antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P15391

Hazard

Information

This CD19 antibody is available for research use only.

Dilution Range

ELISA (order BSA/sodium azide-free format for coating)

Expiration Date

12 months from date of receipt.