

Product Datasheet

EPS8L2 antibody (Biotin) (orb38575)

Description	EPS8L2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Epidermal growth factor receptor kinase substrate 8-like protein 2 protein (201-500AA)
Target	EPS8L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H6S3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.