

Product Datasheet

TRP1 / Tyrosinase-Related Protein-1 Antibody (orb385760)

Description

It reacts with a 75kDa melanocyte-

	Tyrosinase-related protein-1 ...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FACS, IF, IHC-Fr
Immunogen	SK-MEL-23 cells were used as the immunogen for this TRP1 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the TRP1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	The optimal dilution of the TRP1 antibody for each application should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P17643
Hazard Information	This TRP1 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (Frozen): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.