

Product Datasheet

C4d / Complement 4d Antibody (orb385789)

Description

This antibody is specific to Complement 4d (C4d)


Species/Host

bou...

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of Human
ki...

Immunogen

Human protein was used as the immunogen for this C4d antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the C4d antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the C4d antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This C4d antibody is available for research use only.

Entrez

[720](#), [721](#)

Dilution Range

Immunohistochemistry (FFPE): 1:200-1:400 for 30 min at RT

Expiration Date

12 months from date of receipt.