

Product Datasheet

Secretory Component p85 Antibody / ECM1 (orb385823)

Description

Extracellular matrix protein 1 (ECM1), also



protein that i...

Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

IHC-P

Applications


Immunohistochemical
staining of human
Co...

Immunogen

Recombinant human protein was used as the immunogen for the Secretory Component Glycoprotein antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Secretory Component Glycoprotein antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Secretory Component Glycoprotein antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

Q16610

Hazard
Information

This Secretory Component Glycoprotein antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

Expiration Date

12 months from date of receipt.