

Product Datasheet

EMI1 Antibody (Early Mitotic Inhibitor-1) (orb385841)

Description

Early Mitotic Inhibitor-1 regulates mitosis by

complex/cyclosome (...)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested
Applications

IHC-P, WB

Immunogen

A human partial recombinant N-terminal protein fragment was used as the immunogen for the EMI-1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the EMI-1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the EMI-1 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q9UKT4

Hazard
Information


This EMI-1 antibody is available for research use only.

Dilution Range


Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

Expiration Date

12 months from date of receipt.



Western blot analysis
of human
HeLa(lane...



Immunohistochemical
staining of Human
OV...