

Product Datasheet

MUC6 Antibody / Gastric Mucin (orb385877)

Description

Mucin 6 may provide a mechanism for modulation


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

IHC-P

Applications


Immunohistochemical
staining of Human
ga...

Immunogen

Recombinant human protein was used as the immunogen for the Mucin-6 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Mucin-6 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titering of the Mucin-6 antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q6W4X9

Hazard
Information

This Mucin-6 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 2-4ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

Expiration Date

12 months from date of receipt.