

Product Datasheet

Fascin Antibody (orb385897)

Description

Recognizes a protein of 55kDa, which is identified as



Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

FACS, IF, IHC-P, WB

Immunohistochemical
staining of human
Ho...

Immunogen

Full length recombinant human protein was used as the immunogen for the Fascin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage


Store the Fascin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Fascin antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.



Immunohistochemical
staining of Human
OV...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q16658

Hazard

Information


This Fascin antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

Expiration Date

12 months from date of receipt.



Western blot analysis
of HeLa cell lysat...