

Product Datasheet

GRP94 Antibody / HSP90B1 (orb385913)

Description

Heat shock protein 90kDa beta member 1, known


Species/Host

Rat

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of Human
br...

Immunogen

Recombinant full-length human HSP90B1 protein was used as the immunogen for the GRP94 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the GRP94 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titering of the GRP94 antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rat IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
P14625
Hazard
Information

This GRP94 antibody is available for research use only.

Dilution Range

Immunohistology (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

Expiration Date

12 months from date of receipt.