

Product Datasheet


VEGFA Antibody / Vascular Endothelial Growth Factor A (orb385915)

Description

Vascular Endothelial Growth Factor is a



family. It is a hepari...



Immunohistochemical
staining of Human
OV...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Human recombinant protein was used as the immunogen for the VEGF antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the VEGF antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titration of the VEGF antibody may be required for optimal performance.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 8, for 10-20 min followed by cooling at RT for 20 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P15692

Hazard
Information

This VEGF antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 4-8ug/ml for 30 min at RT (1)

Expiration Date

12 months from date of receipt.