

Product Datasheet

HPV-16 L1 Antibody (orb385921)

Description

Infection with specific types of HPV has been



cervic...

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

Applications

IHC-P

Immunohistochemical
staining of Human
ce...

Immunogen

Human papilloma virus type 16, major capsid protein L1, was used as the immunogen for the HPV-16 L1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the HPV-16 L1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titration of the HPV-16 L1 antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P03101

Hazard

Information

This HPV-16 L1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT (1)

Expiration Date

12 months from date of receipt.