


Product Datasheet

LRRC4C antibody (orb38593)


Description	LRRC4C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 4C protein (201-500AA)
Target	LRRC4C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCJ2
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gl...



Immunohistochemical staining of human br...



Immunofluorescence analysis of MCF-7 cel...