

Product Datasheet

Interferon gamma Antibody / IFNG (orb385937)

Description

Recognizes a protein of 20-25kDa, identified as


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical staining of Human resected tissue

Immunogen


Human recombinant protein was used as the immunogen for the IFN gamma antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the IFN gamma antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical staining of Human resected tissue

Note

For research use only

Application notes

Optimal dilution of the IFN gamma antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype


Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography



Immunohistochemical staining of Human resected tissue

Uniprot ID

P01579

Hazard
Information

This IFN gamma antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.