

Product Datasheet

CD15 Antibody / SSEA-1 (orb385947)

Description

CD15 plays a role in mediating phagocytosis,



on...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P

Applications

Immunohistochemical
staining of human
Ho...

Immunogen

Purified human neutrophils were used as the immunogen for the CD15 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CD15 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the CD15 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity

Uniprot ID

P22083

Hazard

This CD15 antibody is available for research use only.

Information
Dilution Range

Flow cytometry: 2-4ug/million cells, Immunofluorescence: 2-4ug/ml, Immunohistochemistry (FFPE): 5-10ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.