

Product Datasheet

Gamma Catenin Antibody / Plakoglobin (orb385949)

Description

Plakoglobin, or gamma Catenin, is a common

plaque...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

Applications

FACS, IF, WB

Immunogen

Recombinant full length chicken gamma Catenin was used as the immunogen for the Plakoglobin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Plakoglobin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting amounts. Titration of the Plakoglobin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P14923

Hazard

Information


This Plakoglobin antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of human HT20 (lane...