

Product Datasheet

**Adipophilin Antibody / Perilipin 2 / ADRP / PLIN2
(orb385951)**

Description

Adipophilin, also called Perilipin 2, ADRP (Adipose


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunohistochemical
staining of Human
ad...

Immunogen

A human partial recombinant protein was used as the immunogen for this ADRP antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

The ADRP antibody (with azide) can be stored at 2-8°C. The azide-free format should be aliquoted and stored at -20°C or colder.

Note

For research use only

Application notes

Titration of the ADRP antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard
Information

This ADRP antibody is available for research use only.

Entrez

[123](#)

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Immunohistochemistry (FFPE): 2-4ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.