


Product Datasheet

Beta Catenin Antibody (orb385953)

Description

The protein encoded by this gene is part of a complex of



Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FACS, IF, WB

Applications

Immunogen

A recombinant fusion protein consisting of the N-terminal half of human beta-Catenin fused to maltose binding protein was used as the immunogen for the Beta Catenin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Beta Catenin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titration of the Beta Catenin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[P35222](#)

Hazard

Information

This Beta Catenin antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.