

Product Datasheet


Alpha-1-Antichymotrypsin Antibody / AACT / SERPINA3 (orb385969)

Description

Alpha 1-antichymotrypsin, or AACT, is a



liver as a sing...



Immunohistochemical
staining of Human
pa...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A human partial recombinant protein was used as the immunogen for this Alpha-1-Antichymotrypsin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Alpha-1-Antichymotrypsin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting amounts. Titration of the Alpha-1-Antichymotrypsin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Protein G affinity

Uniprot ID

P01011

Hazard

Information

This Alpha-1-Antichymotrypsin antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.