

Product Datasheet

Caldesmon / CAD Antibody HMW (orb385973)

Description

Caldesmon is a developmentally regulated



muscle contractio...

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical
staining of Human
ut...

Immunogen

A partial human protein was used as the immunogen for the Caldesmon antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Caldesmon antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical
staining of Rat uter...

Note

For research use only

Application notes

Titering of the Caldesmon antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity chromatography

Uniprot ID

[Q05682](#)

Hazard

Information

This Caldesmon antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.