

Product Datasheet

Gamma Catenin Antibody / Plakoglobin (orb385977)

Description

Plakoglobin, also known as junction plakoglobin



is enc...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical staining of Human ce...

Immunogen


Recombinant human protein was used as the immunogen for the Plakoglobin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Plakoglobin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical staining of Human sk...

Note

For research use only

Application notes

The stated application concentrations are suggested starting amounts. Titration of the Plakoglobin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P14923

Hazard
Information

This Plakoglobin antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.1-0.2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.