

Product Datasheet

AMACR / p504S Antibody (Prostate Cancer Marker) (orb385979)

Description

AMACR, also known as p504S is an enzyme that



oxidatio...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

IHC-P, WB

Immunohistochemical
staining of Human
pr...

Immunogen

A recombinant human protein was used as the immunogen for the AMACR antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the AMACR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the AMACR antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, lambda

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q9UHK6

Hazard

Information

This AMACR antibody is available for research use only.

Dilution Range

Western blot: 0.5-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.