

Product Datasheet

Beta Catenin Antibody (orb385983)

Description

Beta-Catenin associates with the cytoplasmic

funct...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A partial human recombinant protein was used as the immunogen for the Beta Catenin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Beta Catenin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Beta Catenin antibody to be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P35222

Hazard

Information

This Beta Catenin antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 0.1-0.2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of Human to...



Immunohistochemical staining of Human pa...



Immunohistochemical staining of Mouse co...