

Product Datasheet

Beta Catenin Antibody (orb385990)

Description

Beta Catenin is a dual function protein, regulating the



Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

Applications

FACS, IF, IHC-P, WB

Immunohistochemical staining of Human to...

Immunogen

A C-terminal 100 amino acid segment from chicken beta-Catenin was used as the immunogen for the Beta Catenin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Beta Catenin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Beta Catenin antibody to be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P35222

Hazard

Information


This Beta Catenin antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 0.5-1ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Western blot analysis of human A431(lane...