

Product Datasheet

Desmoglein 1 Antibody / DSG1 (orb386001)

Description

Desmoglein 1 (DSG1) is a member of the



intercellular adhe...

Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

IHC-P

Applications


Immunohistochemical
staining of Human
sk...

Immunogen

A recombinant protein corresponding to the intracellular domain of human Desmoglein 1 was used as the immunogen for the Desmoglein 1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Desmoglein 1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Desmoglein 1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q02413](#)

Hazard
Information

This Desmoglein 1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.