

Product Datasheet

Chromogranin A Antibody (orb386014)

Description Chromogranin A is present in neuroendocrine cells


...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

 Immunohistochemical
staining of Human
pa...

Immunogen


Recombinant human protein was used as the immunogen for the Chromogranin A antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Chromogranin A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).


 Immunohistochemical
staining of Human
pa...

Note

For research use only

Application notes

Optimal dilution of the Chromogranin A antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity chromatography

Uniprot ID
P10645
Hazard
Information

This Chromogranin A antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.25-0.5ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.