

Product Datasheet

GCDFP-15 Antibody / PIP / Prolactin-Induced Protein (orb386022)

Description

It recognizes a protein of 15kDa, identified as


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested
Applications

FACS, IF, IHC-P, WB

Immunohistochemical
staining of Human
br...

Immunogen


A human partial recombinant protein corresponding to amino acids 41-146 was used as the immunogen for the PIP antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the PIP antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Western blot analysis
of human pancreas(...

Note

For research use only

Application notes

Titration of the PIP antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P12273

Hazard
Information

This PIP antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml, Western blot: 0.5-1ug/ml, Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.