

Product Datasheet

Galectin 13 Antibody / LGALS13 (orb386036)

Description

It recognizes a 32kDa protein, which is identified as



Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P, WB

Applications

Immunohistochemical staining of Human pl...

Immunogen


A human partial recombinant protein corresponding to amino acids 23-134 was used as the immunogen for the Galectin 13 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Galectin 13 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical staining of Human sp...

Note

For research use only

Application notes

Titering of the Galectin 13 antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q9UHV8](#)

Hazard

Information


This Galectin 13 antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 0.5-1ug/ml, Western blot: 0.5-1ug/ml, Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Western blot analysis of human partial r...