

Product Datasheet

CD80 Antibody (orb386038)

Description

T cell proliferation and lymphokine production

	antigen, fo...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant human protein used as the immunogen for the CD80 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the CD80 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the CD80 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P33681
Hazard Information	This CD80 antibody is available for research use only.
Dilution Range	ELISA (order BSA/sodium azide-free format for coating)
Expiration Date	12 months from date of receipt.